

BE BEAUTIFUL.  
BE EXTRAORDINARY.  
BE MIXO.

+ (82) 1644-2831

55, Convensia-daero, of. 202, Yeonsu-gu,

Incheon, South Korea

DERMAKOR

BY DERMAKOR

# MIXO

Sterile Absorbable  
Hyaluronic Acid Dermal Filler



BE BEAUTIFUL.

BE EXTRAORDINARY.

BE MIXO.

BE BEAUTIFUL.

BE EXTRAORDINARY.

BE MIXO.

## Why MIXO?

**M** means Beauty from Korean Hanja and **XO** is from eXtraOrdinary.

**MIXO** is premium monophasic hyaluronic acid filler for beautiful and unique people who want to emphasize their natural beauty and individuality

What is Hyaluronic Acid Filler?

Safety, Biodegradable and biocompatible material

HA has a much greater capacity to retain than any other natural or synthetic polymer. Hydrates, regenerates skin, improves elastin and collagen production...

Reduces the appearance of fine wrinkles and returns lost contours of the face.

Can be removed easily and perfectly.



**ENHANCED NATURAL BEAUTY AND PERSONALITY WITH MIXO**

• **MIXO** Produced by patented Nano Mesh Cross-Linking Technology for maintaining safety and long lasting filler.



Raw material HA, tangled and uneven texture



Even and stable structure, ready for smooth cross-linking



Cross-linked HA, same with other dermal fillers



Strengthened and upgraded structure of cross-linked MIXO

“ **MIXO** Innovative Technology produces extremely small and solid molecules, which makes more comfortable and painless treatment. **MIXO** and fluid nature of HA allows easy penetration into the skin maximizing effect of reduced wrinkles and higher satisfaction from the procedure. ”

• **Other fillers** EDCE Cross-linking to Big, Middle, Small particles and mixing Cross-linked Gel with solution.



Cross-link



cross linked HA Mashing to the micro-size



HA in Water



Finished Product

**MIXO**

DVS makes Micro-Beads to molecular size without mixing with solution.



Cross-link



Cross-linking



Micro-Gel Aggregation



Final Product

## HIGH VISCOELASTICITY AND HIGH COHESIVENESS

- Patented manufacturing technology make possible to form a MASS GEL that best represents the performance of the supporting structure in the human body, High viscosity provides necessary volume and beautiful shapes



	MIXO	Product A	Product B	Product
Viscosity	800	380	610	520

- High cross-linking rate possess long-lasting effect.

MIXO	A	B	C




## SAFETY AND BIODEGRADABLE

### - Refractive absorbance analysis

raw material




### - MIXO Complete product



\* The wavelength pattern of MIXO is the same as its raw material showing that the product has not altered its physical properties.

- MIXO is a Sterile Absorbable Medical Device that can be safely and perfectly dissolved through Hyaluronidase.



## Part of Procedure

- MIXO LIGHT PLUS
- MIXO DEEP PLUS
- MIXO VOLUME PLUS



## PRODUCT INFORMATION

	MIXO LIGHT	MIXO DEEP	MIXO VOLUME	MIXO LIGHT Plus	MIXO DEEP Plus	MIXO VOLUME Plus
COMPOSITION	HA 20mg/mL	HA 20mg/mL	HA 20mg/mL	HA 20mg/mL, Lidocaine 0.2%	HA 20mg/mL, Lidocaine 0.2%	HA 20mg/mL, Lidocaine 0.2%
Packaging Unit	1.2mL x 1	1.0mL x 1	1.2mL x 1	1.1mL x 1	1.1mL x 1	1.1mL x 1
Needle Type	30G x 2pcs	27G x 2pcs	25G x 2pcs	30G x 2pcs	27G x 2pcs	25G x 2pcs
Storage	2°C-25°C	2°C-25°C	2°C-25°C	2°C-25°C	2°C-25°C	2°C-25°C

“ Product name: Sterile absorbable Hyaluronic Acid Dermal Filler ”

Intended Use: Temporary improvement of the facial wrinkles

Sterile Medical Device/ Single Use Only/ Do not reuse

Storage Condition: Store at room temperature, Avoid direct sunlight,

Please refer to instruction for use for more details.

